

Anti-Serpinf2 Rabbit Polyclonal Antibody For Wb And Elisa - P08697

Item Number: CM0013045

Introduction

High-quality rabbit polyclonal antibody against SERPINF2 (alpha-2-antiplasmin). Validated for WB and ELISA, cross-reacts with human, mouse, rat. Ideal for fibrinolysis and serpin research. - P08697

[Learn More](#)

Parameter	Value
Product Name	SERPINF2 Rabbit pAb
Remarks/Alias	AAP, API, PLI, A2AP, alpha2AP, ALPHA-2-PI, SERPINF2
Species	Human
Gene ID (Human)	5345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 302-491 of human SERPINF2 (NP_001159392.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08697
Protein Weight	55kDa
Shipping	Ice bag