


Anti-Serping1 Rabbit Monoclonal Antibody For Wb And Elisa - P05155

Item Number: CM0002798



Introduction

Rabbit monoclonal antibody against human SERPING1 (Plasma protease C1 inhibitor), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for complement pathway studies.

[Learn More](#)

Parameter	Value
Product Name	SERPING1 Rabbit mAb
Remarks/Alias	C1IN; C1NH; HAE1; HAE2; C1INH; SERPING1
Species	Human
Gene ID (Human)	710
Gene ID	710
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human SERPING1 (P05155).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05155
Protein Weight	55kDa
Shipping	Ice bag