


Anti-Serping1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05155

Item Number: CM0016494



Introduction

Rabbit polyclonal antibody targeting human SERPING1 (C1 inhibitor). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for complement pathway, coagulation, and kinin system research.

[Learn More](#)

Parameter	Value
Product Name	SERPING1 Rabbit pAb
Remarks/Alias	C1IN, C1NH, HAE1, HAE2, C1INH, SERPING1
Species	Human
GeneID	710
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-500 of human SERPING1 (NP_000053.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05155
Protein Weight	55kDa
Shipping	Ice bag