

Anti-Serpini1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99574

Item Number: CM0013327



Introduction

Rabbit polyclonal antibody targeting human SERPINI1/neuroserpin, suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying synaptic plasticity, neurodegeneration, and protease inhibition in brain tissues.

[Learn More](#)

Parameter	Specification
Product Name	SERPINI1 Rabbit pAb
Remarks/Alias	PI12; HNS-S1; HNS-S2; neuroserpin; SERPINI1
Species	Human
GeneID (human)	5274
GeneID	5274
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 181-410 of human SERPINI1 (NP_005016.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q99574
Protein Weight	46kDa
Shipping	Ice bag