

Anti-Serpini2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - 075830

Item Number: CM0013626

Introduction

SERPINI2 rabbit polyclonal antibody validated for Western blot, immunofluorescence (IF-P), and ELISA. Detects human SERPINI2 and shows cross-reactivity with mouse and rat, supporting protease inhibitor and pancreatic research.

[Learn More](#)

Parameter	Value
Product Name	SERPINI2 Rabbit pAb
Remarks/Alias	MEPI; PI14; PANCPI; TSA2004; SERPINI2
Species	Human
GeneID (Human)	5276
GeneID	5276
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 56-405 of human SERPINI2 (NP_006208.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	075830
Protein Weight	46kDa
Shipping	Ice bag