


## Anti-Sesn2 Rabbit Polyclonal Antibody For Wb, Elisa - P58004

Item Number: CM0016378



### Introduction

Rabbit polyclonal antibody targeting human SESN2 (Sestrin-2/HI95), validated for WB and ELISA. Cross-reacts with mouse. Ideal for mTOR signaling and oxidative stress research.

[Learn More](#)

| Parameter        | Value                    |
|------------------|--------------------------|
| Product Name     | SESN2 Rabbit pAb         |
| Remarks/Alias    | HI95; SES2; SEST2; SESN2 |
| Species          | Human                    |
| Gene ID (Human)  | 83667                    |
| Immunogen        | Recombinant protein      |
| Source           | Rabbit                   |
| Category         | Polyclonal Antibodies    |
| Application      | WB, ELISA                |
| Cross Reactivity | Human, Mouse             |
| SWISS            | P58004                   |
| Protein Weight   | 54kDa                    |
| Shipping         | Ice bag                  |