

Anti-Sesn3 Polyclonal Antibody For Wb, Elisa - P58005

Item Number: CM0025170

Introduction

Rabbit polyclonal anti-SESN3 antibody validated for WB and ELISA. Detects human SESN3 (57kDa), a leucine sensor and oxidative stress protector. Ideal for metabolic and stress research. - P58005

[Learn More](#)

Specification	Details
Product Name	SESN3 Rabbit pAb
Remarks/Alias	SEST3; SESN3
Species	Human
Gene ID (Human)	143686
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human SESN3 (NP_653266.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P58005
Protein Weight	57kDa
Shipping	Ice bag