

Anti-Sestrin 1 (Pa26/Sest1) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Y6P5

Item Number: CM0002237

Introduction

Rabbit monoclonal antibody against human sestrin 1 (PA26/SEST1), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse. Ideal for studying TORC1 signaling, oxidative stress, and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	sestrin 1 Rabbit mAb
Remarks/Alias	PA26; SEST1; sestrin 1
Species	Human
Gene ID	27244
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human sestrin 1 (Q9Y6P5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6P5
Protein Weight	57kDa
Shipping	Ice bag