

Anti-Sestrin 2 Rabbit Monoclonal Antibody For Wb, Elisa - P58004

Item Number: CM0007304



Introduction

Rabbit monoclonal Sestrin 2 antibody validated for WB and ELISA, cross-reacting with human, mouse, and rat SESN2. This stress-inducible metabolic regulator is a leucine sensor modulating mTORC1; suitable for oxidative stress and nutrient signaling research.

[Learn More](#)

Parameter	Value
Product Name	Sestrin 2 Rabbit mAb
Remarks/Alias	HI95; SES2; SEST2; Sestrin 2
Species	Human
GeneID (human)	83667
Immunogen	Recombinant fusion protein containing amino acids 300-400 of human Sestrin 2 (NP_113647.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P58004
Protein Weight	54kDa
Shipping	Ice bag