


Anti-Set1A Polyclonal Antibody For Wb, Elisa - O15047

Item Number: CM0022164



Introduction

High-quality rabbit polyclonal antibody targeting human SET1A (O15047), validated for WB and ELISA; cross-reacts with human and mouse. Ideal for chromatin remodeling and histone methylation research.

[Learn More](#)

Core Parameters	Details
Product Name	SET1A Rabbit pAb
Remarks/Alias	Set1; EPEDD; KMT2F; Set1A; NEDSID; SET1A
Species	Human
GeneID (Human)	9739
GeneID	9739
Immunogen	Recombinant fusion protein containing amino acids 258-458 of human SET1A (O15047)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15047
Protein Weight	186kDa
Shipping	Ice bag