

Anti-Setbp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6X0

Item Number: CM0030226

Introduction

High-specificity anti-SETBP1 rabbit polyclonal antibody for Western blot, immunofluorescence, and ELISA. Reacts with human SETBP1 (175kDa). Useful for studying myeloid leukemias and intellectual disability.

[Learn More](#)

Parameter	Value
Product Name	SETBP1 Rabbit pAb
Remarks/Alias	SEB; MRD29; SETBP1
Species	Human
Gene ID (Human)	26040
Gene ID	26040
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-242 of human SETBP1 (NP_001123582.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y6X0
Protein Weight	175kDa
Shipping	Ice bag