

Anti-Setd1B Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ups6

Item Number: CM0019895



Introduction

Rabbit polyclonal antibody targeting human SETD1B (KMT2G), suitable for Western blotting (WB) and ELISA. Reacts with human, mouse, and rat samples. For chromatin biology and hematopoiesis research.

[Learn More](#)

Parameter	Value
Product Name	SETD1B Rabbit pAb
Remarks/Alias	KMT2G; Set1B; IDDSELD; SETD1B
Species	Human
GeneID (Human)	23067
GeneID	23067
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1673-1770 of human SETD1B (NP_001340274.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPS6
Protein Weight	213kDa
Shipping	Ice bag