


Anti-Setd2 Polyclonal Antibody For Wb, Elisa - Q9Byw2

Item Number: CM0018249



Introduction

Rabbit polyclonal antibody targeting human SETD2, a histone methyltransferase. Validated for WB and ELISA, with cross-reactivity to human. Suitable for epigenetic, cancer, and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	SETD2 Rabbit pAb
Remarks/Alias	LLS, HYPB, SET2, HIF-1, HIP-1, KMT3A, MRD70, RAPAS, HBP231, HSPC069, p231HBP, SETD2
Species	Human
GeneID (Human)	29072
GeneID	29072
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 803-1103 of human SETD2 (NP_054878.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BYW2
Protein Weight	288kDa
Shipping	Ice bag