


Anti-Setd2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q9Byw2

Item Number: CM0029520



Introduction

High-quality anti-SETD2 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Targets histone methyltransferase involved in H3K36me3, DNA repair, and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	SETD2 Rabbit pAb
Remarks/Alias	LLS; HYPB; SET2; HIF-1; HIP-1; KMT3A; MRD70; RAPAS; HBP231; HSPC069; p231HBP; SETD2
Species	Human
GeneID (Human)	29072
GeneID	29072
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 2475-2564 of human SETD2 (NP_054878.5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYW2
Protein Weight	288kDa
Shipping	Ice bag