

Anti-Setd3 Polyclonal Antibody For Wb, Ip, Elisa - Q86Tu7

Item Number: CM0018943



Introduction

Rabbit polyclonal antibody against human SETD3 (hSETD3/C14orf154), validated for WB, IP, ELISA. Specific for actin-histidine N-methyltransferase, ~67 kDa. Ideal for cytoskeletal and methylation studies.

[Learn More](#)

Parameter	Value
Product Name	SETD3 Rabbit pAb
Remarks/Alias	hSETD3; C14orf154; SETD3
Species	Human
GeneID (Human)	84193
GeneID	84193
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 72-329 of human SETD3 (NP_115609.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q86TU7
Protein Weight	67kDa
Shipping	Ice bag