

Anti-Setd5 Rabbit Polyclonal Antibody For Wb - Q9C0A6

Item Number: CM0026916

Introduction

Rabbit polyclonal antibody targeting human SETD5, validated for WB, IHC-P, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Ideal for studying chromatin regulation in brain development.

[Learn More](#)

Parameter	Value
Product Name	SETD5 Rabbit pAb
Remarks/Alias	MRD23; SETD5A; SETD5
Species	Human
GeneID (Human)	55209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1054-1153 of human SETD5 (NP_001073986.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C0A6
Protein Weight	158kDa
Shipping	Ice bag