

Anti-Setd6 Polyclonal Antibody For Wb, Elisa - Q8Tbk2

Item Number: CM0027756



Introduction

Rabbit polyclonal antibody against human SETD6, an N-lysine methyltransferase involved in NF-κB signaling and stem cell self-renewal. Validated for WB and ELISA; cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	SETD6 Rabbit pAb
Remarks/Alias	SETD6
Species	Human
GeneID (Human)	79918
GeneID	79918
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 301-473 of human SETD6 (NP_001153777.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TBK2
Protein Weight	53kDa
Shipping	Ice bag