


Anti-Setdb1 Rabbit Polyclonal Antibody (Ko Validated) For Wb, If-P, Elisa - Q15047

Item Number: CM0028862



Introduction

Rabbit polyclonal antibody against histone methyltransferase SETDB1, KO-validated. Suitable for WB, IF-P, ELISA in human, mouse, rat. Ideal for studying transcriptional repression, heterochromatin formation, and epigenetic mechanisms in cancer and stem cells.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SETDB1 Rabbit pAb
Remarks/Alias	ESET; KG1T; KMT1E; TDRD21; H3-K9-HMTase4; B1
Species	Human
GeneID (Human)	9869
GeneID	9869
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-397 of human SETDB1 (NP_001230420.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15047
Protein Weight	143kDa
Shipping	Ice bag