

Anti-Setdb2 Rabbit Polyclonal Antibody For Wb, Elisa - Q96T68

Item Number: CM0029905



Introduction

High-quality rabbit polyclonal antibody targeting human SETDB2 (KMT1F), validated for Western blot and ELISA. Cross-reacts with mouse and rat, making it ideal for epigenetic research on histone methylation and gene silencing.

[Learn More](#)

Parameter	Value
Product Name	SETDB2 Rabbit pAb
Remarks/Alias	CLLD8; CLLL8; KMT1F; C13orf4; SETDB2
Species	Human
GeneID (Human)	83852
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 148-378 of human SETDB2 (NP_114121.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96T68
Protein Weight	82kDa
Shipping	Ice bag