

Anti-Setmar Rabbit Polyclonal Antibody For IHC-P And Elisa - Q53H47

Item Number: CM0025396



Introduction

Rabbit polyclonal antibody against human SETMAR/Metnase for IHC-P and ELISA, targeting a 78kDa histone methyltransferase-transposase fusion involved in DNA repair.

[Learn More](#)

Parameter	Value
Product Name	SETMAR Rabbit pAb
Remarks/Alias	Mar1; METNASE; SETMAR
Species	Human
Gene ID	6419
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human SETMAR (NP_006506.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q53H47
Protein Weight	78kDa
Shipping	Ice bag