

Anti-Set/Taf1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q01105

Item Number: CM0008353



Introduction

Rabbit monoclonal antibody targeting human SET/TAF1 protein. Validated for WB, IHC-P, IF/ICC, IP, and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for apoptosis, transcription, and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	SET/TAF1 Rabbit mAb
Remarks/Alias	2PP2A; IGAAD; MRD58; TAF-I; I2PP2A; IPP2A2; PHAPII; TAF-IBETA; SET/TAF1
Species	Human
GeneID (human)	6418
GeneID	6418
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-290 of human SET/TAF1 (Q01105).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01105
Protein Weight	33kDa
Shipping	Ice bag