


Anti-Sf3A1 Rabbit Polyclonal Antibody For Wb, Elisa - Q15459

Item Number: CM0028822



Introduction

High-quality anti-SF3A1 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat. Targets splicing factor 3A subunit 1, essential for spliceosome assembly and pre-mRNA splicing.

[Learn More](#)

Parameter	Value
Product Name	SF3A1 Rabbit pAb
Remarks/Alias	PRP21; PRPF21; SAP114; SF3A120; SF3A1
Species	Human
GeneID (human)	10291
GeneID	10291
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 715-785 of human SF3A1 (Q15459)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q15459
Protein Weight	89kDa
Shipping	Ice bag