


## Anti-Sf3A3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q12874

Item Number: CM0028032



### Introduction

Anti-SF3A3 polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat SF3A3, a 59kDa splicing factor subunit of U2 snRNP. Suitable for mRNA processing research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | SF3A3 Rabbit pAb                                                                                                   |
| Remarks/Alias    | PRP9; PRPF9; SAP61; SF3a60; SF3A3                                                                                  |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 10946                                                                                                              |
| Gene ID          | 10946                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 292-501 of human SF3A3 (NP_006793.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q12874                                                                                                             |
| Protein Weight   | 59kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |