

Anti-Sf3B1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - 075533

Item Number: CM0011698



Introduction

Rabbit monoclonal antibody targeting human SF3B1, the U2 snRNP splicing factor. Validated for WB, IHC-P, IF/ICC, and ELISA, with cross-reactivity to mouse and rat. Suitable for spliceosome and pre-mRNA processing research.

[Learn More](#)

Parameter	Value
Product Name	SF3B1 Rabbit mAb
Remarks/Alias	MDS; PRP10; Hsh155; PRPF10; SAP155; SF3b155; SF3B1
Species	Human
Gene ID (Human)	23451
Gene ID	23451
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SF3B1 (075533).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075533
Protein Weight	146kDa
Shipping	Ice bag