

Anti-Sf3B1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 075533

Item Number: CM0023426



Introduction

SF3B1 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA; detects human, mouse, rat SF3B1 (splicing factor 3B subunit 1, ~146kDa). Ideal for spliceosome and pre-mRNA splicing research.

[Learn More](#)

Parameter	Value
Product Name	SF3B1 Rabbit pAb
Remarks/Alias	MDS; PRP10; Hsh155; PRPF10; SAP155; SF3b155; SF3B1
Species	Human
GeneID (Human)	23451
GeneID	23451
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human SF3B1 (NP_001005526.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075533
Protein Weight	146kDa
Shipping	Ice bag