

Anti-Sf3B3/Sap130 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15393

Item Number: CM0008427



Introduction

Rabbit anti-human SF3B3/SAP130 monoclonal antibody validated in WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat SF3B3. Suitable for splicing research and spliceosome studies.

[Learn More](#)

Parameter	Value
Product Name	SF3B3/SAP130 Rabbit mAb
Remarks/Alias	RSE1; SAP130; SF3b130; STAF130; SF3B3/SAP130
Species	Human
GeneID (Human)	23450
GeneID	23450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1101-1217 of human SF3B3/SAP130 (Q15393)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15393
Protein Weight	136kDa
Shipping	Ice bag