

Anti-Sf3B6 Polyclonal Antibody For Wb And Elisa - Q9Y3B4

Item Number: CM0018042



Introduction

SF3B6 Rabbit pAb validated for WB and ELISA, detecting human and mouse SF3B6 (14 kDa), a core U2 snRNP component essential for spliceosome assembly and branch site recognition.

[Learn More](#)

Parameter	Value
Product Name	SF3B6 Rabbit pAb
Remarks/Alias	P14; Ht006; SAP14; SAP14a; SF3B14; CGI-110; HSPC175; SF3B14a; SF3B6
Species	Human
Gene ID (Human)	51639
Gene ID	51639
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human SF3B6 (NP_057131.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y3B4
Protein Weight	15kDa
Shipping	Ice bag