


Anti-Sfpq Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P23246

Item Number: CM0004430



Introduction

High-quality anti-SFPQ rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Targets splicing factor proline- and glutamine-rich involved in pre-mRNA splicing, transcriptional regulation, and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	SFPQ Rabbit mAb
Remarks/Alias	PSF; POMP100; PPP1R140; SFPQ
Species	Human
Gene ID	6421
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 600-707 of human SFPQ (P23246).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23246
Protein Weight	76kDa
Shipping	Ice bag