


Anti-Sfpq Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P23246

Item Number: CM0013886



Introduction

Rabbit polyclonal antibody against human SFPQ (PSF), validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Molecular weight: 76 kDa. Suitable for splicing and transcription research.

[Learn More](#)

Parameter	Value
Product Name	SFPQ Rabbit pAb
Remarks/Alias	PSF; POMP100; PPP1R140; SFPQ
Species	Human
Gene ID (Human)	6421
Gene ID	6421
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 578-707 of human SFPQ (NP_005057.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23246
Protein Weight	76kDa
Shipping	Ice bag