

Anti-Sfrp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N474

Item Number: CM0025081

Introduction

SFRP1 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in human and mouse. Targets secreted frizzled-related protein 1, a Wnt signaling modulator, for cell growth and differentiation research.

[Learn More](#)

Parameter	Value
Product Name	SFRP1 Rabbit pAb
Remarks/Alias	FRP; FRP1; FrzA; FRP-1; SARP2; SFRP1
Species	Human
Gene ID	6422 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-314 of human SFRP1 (NP_003003.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N474
Protein Weight	35kDa
Shipping	Ice bag