

Anti-Sfrp2 Polyclonal Antibody For Wb And Elisa - Q96Hf1

Item Number: CM0027405

Introduction

Rabbit polyclonal antibody against human SFRP2, validated for Western blot and ELISA applications. Cross-reacts with mouse; ideal for Wnt signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	SFRP2 Rabbit pAb
Remarks/Alias	FRP-2; SARP1; SDF-5; SFRP2
Species	Human
GeneID	6423
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-295 of human SFRP2 (NP_003004.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96HF1
Protein Weight	33kDa
Shipping	Ice bag