

## Anti-Sfrp4 Monoclonal Antibody For Wb, If/Icc, Elisa - Q6Fhj7

Item Number: CM0004919

### Introduction

Recombinant rabbit monoclonal antibody targeting human SFRP4 (sFRP-4), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for Wnt signaling and cancer studies.

[Learn More](#)

| Parameter        | Value                                                                                                          |
|------------------|----------------------------------------------------------------------------------------------------------------|
| Product Name     | SFRP4 Rabbit mAb                                                                                               |
| Remarks/Alias    | PYL; FRP-4; FRPHE; FRZB-2; sFRP-4; SFRP4                                                                       |
| Species          | Human                                                                                                          |
| GeneID (Human)   | 6424                                                                                                           |
| GeneID           | 6424                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 201-300 of human SFRP4 (Q6FHJ7). |
| Source           | Rabbit                                                                                                         |
| Category         | Monoclonal Antibodies                                                                                          |
| Application      | WB, IF/ICC, ELISA                                                                                              |
| Cross Reactivity | Human, Mouse, Rat                                                                                              |
| SWISS            | Q6FHJ7                                                                                                         |
| Protein Weight   | 40 kDa                                                                                                         |
| Shipping         | Ice bag                                                                                                        |