

Anti-Sfrp4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Fhj7

Item Number: CM0026254

Introduction

SFRP4 Rabbit Polyclonal Antibody for WB, IF/ICC and ELISA. Confirmed reactivity with Human, Mouse, Rat. A valuable tool for investigating Wnt pathway modulation and bone development.

[Learn More](#)

Parameter	Value
Product Name	SFRP4 Rabbit pAb
Remarks/Alias	PYL; FRP-4; FRPHE; FRZB-2; sFRP-4; SFRP4
Species	Human
GeneID (Human)	6424
GeneID	6424
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-346 of human SFRP4 (NP_003005.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6FHJ7
Protein Weight	40kDa
Shipping	Ice bag