

Anti-Sfrp5 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q5T4F7

Item Number: CM0021966

Introduction

Rabbit polyclonal antibody against human SFRP5 (SARP3, FRP-1b). Validated for WB, IF-P, IHC-P, ELISA and cross-reacts with human, mouse, rat. Ideal for Wnt signaling and developmental studies.

[Learn More](#)

Parameter	Details
Product Name	SFRP5 Rabbit pAb
Remarks/Alias	SARP3; SFRP5
Species	Human
GeneID (Human)	6425
GeneID	6425
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-317 of human SFRP5 (NP_003006.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5T4F7
Protein Weight	36kDa
Shipping	Ice bag