

Anti-Sfrs9 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13242

Item Number: CM0013054

Introduction

High-quality rabbit polyclonal antibody against human splicing factor SFRS9 (SRp30c). Validated for Western blot (WB), immunofluorescence (IF/ICC), and ELISA. Detects both human and mouse SFRS9. Ideal for alternative splicing and gene expression research.

[Learn More](#)

Parameter	Value
Product Name	SFRS9 Rabbit pAb
Remarks/Alias	SFRS9; SRp30c
Species	Human
Gene ID (Human)	8683
Gene ID	8683
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13242
Protein Weight	26kDa
Shipping	Ice bag