

Anti-Sfswap Rabbit Polyclonal Antibody For Wb And Elisa - Q12872

Item Number: CM0021659



Introduction

Anti-SFSWAP Rabbit Polyclonal Antibody validated for Western blot and ELISA. Targets human splicing factor SFSWAP (SWAP/SFRS8), a regulator of alternative splicing. Recombinant immunogen. For RNA processing research.

[Learn More](#)

Parameter	Value
Product Name	SFSWAP Rabbit pAb
Remarks/Alias	SWAP, SFRS8, SFSWAP
Species	Human
Gene ID (Human)	6433
Gene ID	6433
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-280 of human SFSWAP (NP_004583.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q12872
Protein Weight	105kDa
Shipping	Ice bag