

Anti-Sftpa1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Iwl2

Item Number: CM0029317

Introduction

High-quality anti-SFTPA1 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for pulmonary surfactant research, pathogen recognition studies, and alveolar proteinosis. Reliable tool for innate immunity and lung function.

[Learn More](#)

Parameter	Value
Product Name	SFTPA1 Rabbit pAb
Remarks/Alias	SPA, ILD1, PSAP, PSPA, SP-A, SPA1, PSP-A, SFTP1, SP-A1, COLEC4, SFTPA1B, SP-A1 beta, SP-A1 delta, SP-A1 gamma, SP-A1 epsilon, SFTPA1
Species	Human
GeneID (Human)	653509
GeneID	653509
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 122-221 of human SFTPA1 (NP_005402.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IWL2
Protein Weight	26kDa
Shipping	Ice bag