

Anti-Sftpa1 Rabbit Polyclonal Antibody For Wb - Q8Iwl2/Q8Iwl1

Item Number: CM0024061

Introduction

Rabbit polyclonal antibody against human SFTPA1, validated for Western blot. Cross-reacts with Rat. Ideal for lung surfactant, innate immunity, and respiratory disease research.

[Learn More](#)

Parameter	Value
Product Name	SFTPA1 Rabbit pAb
Remarks/Alias	SPA; ILD1; PSAP; PSPA; SP-A; SPA1; PSP-A; SFTP1; SP-A1; COLEC4; SFTPA1B; SP-A1 beta; SP-A1 delta; SP-A1 gamma; SP-A1 epsilon; SFTPA1
Species	Human
GeneID (Human)	653509/729238
GeneID	653509/729238
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 100-248 of human SFTPA1 (NP_001158116.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Rat
SWISS	Q8IWL2/Q8IWL1
Protein Weight	26kDa
Shipping	Ice bag