


Anti-Sftpc Rabbit Polyclonal Antibody For Wb, Ihc, If, Elisa - P11686

Item Number: CM0020357



Introduction

High-quality rabbit polyclonal antibody against human SFTPC (surfactant protein C) for WB, IHC, IF, and ELISA. Reacts with human, mouse, rat. 21 kDa. Ideal for pulmonary disease research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	SFTPC Rabbit pAb
Remarks/Alias	SP5; SP-C; PSP-C; SFTP2; SMDP2; BRICD6; SFTPC
Species	Human
Gene ID (Human)	6440
Gene ID	6440
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 59-197 of human SFTPC (NP_003009.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11686
Protein Weight	21kDa
Shipping	Ice bag