


Anti-Sftpc Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P11686

Item Number: CM0008894



Introduction

High-specificity rabbit recombinant monoclonal antibody against human SFTPC (surfactant protein C). Validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse. MW 21 kDa.

[Learn More](#)

Parameter	Value
Product Name	SFTPC Rabbit PolymAb®
Remarks/Alias	SP5; SP-C; PSP-C; SFTP2; SMDP2; BRICD6
Species	Human
Gene ID (Human)	6440
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 91-190 of human SFTPC (NP_003009.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11686
Protein Weight	21kDa
Shipping	Ice bag