


Anti-Sftpd (Surfactant Protein D) Rabbit Monoclonal Antibody For Wb, Elisa - P35247

Item Number: CM0007366



Introduction

Rabbit monoclonal antibody targeting human SP-D, validated for Western blot and ELISA. Cross-reacts with mouse and rat. An essential tool for studying pulmonary innate immunity and surfactant biology.

[Learn More](#)

Parameter	Value
Product Name	SFTPD Rabbit mAb
Remarks/Alias	SP-D; PSP-D; SFTP4; COLEC7; SFTPD
Species	Human
GeneID (Human)	6441
GeneID	6441
Immunogen	Synthetic peptide - A synthesized peptide derived from human Surfactant protein D/SP-D
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P35247
Protein Weight	38kDa
Shipping	Ice bag