

Anti-Sfxn1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H9B4

Item Number: CM0014195



Introduction

This rabbit polyclonal antibody targets human SFXN1, a mitochondrial serine transporter crucial for one-carbon metabolism. Validated for WB and ELISA, it cross-reacts with mouse and rat. Ideal for research in amino acid transport and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	SFXN1 Rabbit pAb
Remarks/Alias	TCC; SLC56A1; SFXN1
Species	Human
Gene ID (Human)	94081
Gene ID	94081
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 87-322 of human SFXN1 (NP_073591.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9B4
Protein Weight	36kDa
Shipping	Ice bag