

Anti-Sfxn2 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Nb2

Item Number: CM0020882

Introduction

CamelBio's anti-SFXN2 rabbit polyclonal antibody detects the 36kDa SFXN2 protein in WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for mitochondrial serine transport and iron homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	SFXN2 Rabbit pAb
Remarks/Alias	SLC56A2; SFXN2
Species	Human
GeneID (Human)	118980
GeneID	118980
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-90 of human SFXN2 (NP_849189.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96NB2
Protein Weight	36kDa
Shipping	Ice bag