

Anti-Sfxn3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bwm7

Item Number: CM0017214

Introduction

Rabbit polyclonal antibody targeting human sideroflexin-3 (SFXN3), a mitochondrial serine transporter. Suitable for WB and ELISA applications in one-carbon metabolism and mitochondrial transport studies.

[Learn More](#)

Parameter	Value
Product Name	SFXN3 Rabbit pAb
Remarks/Alias	SFX3; SLC56A3; BA108L7.2; SFXN3
Species	Human
Gene ID (Human)	81855
Gene ID	81855
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human SFXN3 (NP_112233.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BWM7
Protein Weight	36kDa
Shipping	Ice bag