

Anti-Sgcb (Beta-Sarcoglycan) Rabbit Polyclonal Antibody For Wb, Elisa - Q16585

Item Number: CM0029189

Introduction

Rabbit polyclonal antibody against human SGCB (beta-sarcoglycan). Validated for WB and ELISA, it recognizes human, mouse, and rat. Ideal for studying sarcoglycan complex in muscular dystrophy research.

[Learn More](#)

Parameter	Value
Product Name	SGCB Rabbit pAb
Remarks/Alias	A3b; SGC; LGMD2E; LGMDR4; SGCB
Species	Human
GeneID (Human)	6443
GeneID	6443
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 94-318 of human SGCB (NP_000223.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16585
Protein Weight	35kDa
Shipping	Ice bag