

Anti-Sgce Rabbit Polyclonal Antibody For Wb, Elisa - 043556

Item Number: CM0029784

Introduction

High-quality rabbit polyclonal antibody against human SGCE (epsilon-sarcoglycan), validated for Western blot and ELISA. Cross-reacts with mouse. Ideal for studying the dystrophin-glycoprotein complex and its role in muscular dystrophy and neurological disorders.

[Learn More](#)

Product Name	SGCE Rabbit pAb
Remarks/Alias	ESG; DYT11; epsilon-SG; SGCE
Species	Human
GeneID (Human)	8910
GeneID	8910
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-317 of human SGCE (NP_003910.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43556
Protein Weight	50kDa
Shipping	Ice bag