

Anti-Sgk1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00141

Item Number: CM0013455



Introduction

Rabbit polyclonal antibody targeting human SGK1 (serum/glucocorticoid-regulated kinase 1). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for studying ion channel regulation and stress response.

[Learn More](#)

Parameter	Value
Product Name	SGK1 Rabbit pAb
Remarks/Alias	SGK; SGK1
Species	Human
GeneID (human)	6446
GeneID	6446
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 172-431 of human SGK1 (NP_005618.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O00141
Protein Weight	49kDa
Shipping	Ice bag