


Anti-Sgk3 Polyclonal Antibody For Wb, Elisa - Q96Br1

Item Number: CM0026687



Introduction

Rabbit polyclonal antibody against human SGK3 (CISK/SGKL), reactive with human and mouse. Suitable for Western blotting and ELISA. Target is a serine/threonine kinase regulating ion channels and cell growth. 57 kDa band detected.

[Learn More](#)

Parameter	Value
Product Name	SGK3 Rabbit pAb
Remarks/Alias	CISK; SGK2; SGKL; SGK3
Species	Human
GeneID (Human)	23678
GeneID	23678
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human SGK3 (NP_037389.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96BR1
Protein Weight	57kDa
Shipping	Ice bag