

Anti-Sgpl1 Polyclonal Antibody For If-P, Ihc-P, Elisa - O95470

Item Number: CM0026481



Introduction

Rabbit polyclonal antibody against human SGPL1 (S1P lyase). Validated for IF-P, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for sphingolipid metabolism and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	SGPL1 Rabbit pAb
Remarks/Alias	SPL; S1PL; NPHS14
Species	Human
GeneID (human)	8879
GeneID	8879
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95470
Protein Weight	64kDa
Shipping	Ice bag