

Anti-Sgpl1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95470

Item Number: CM0015790



Introduction

Rabbit polyclonal antibody against human SGPL1 (Sphingosine-1-phosphate lyase 1). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for lipid metabolism and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	SGPL1 Rabbit pAb
Remarks/Alias	SPL; S1PL; NPHS14; SGPL1
Species	Human
GeneID (Human)	8879
GeneID	8879
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-390 of human SGPL1 (NP_003892.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95470
Protein Weight	64kDa
Shipping	Ice bag